



## DIGITAL GAS FLOWMETER (VKIT-GFM)



Supplied with 6-Point NIST Traceable Calibration upto 1000ml/min.

Ideal for all GC Flowrates, including column flow, split, purge, fuel & make-up

### **Specifications:**

- Measuring Range: 0.01 – 1000ml/min.
- Display: LCD
- Accuracy:  $\pm 3\%$  of reading, or 0.2ml/min., whichever is greater
- Standard calibration points: 5, 50, 100, 500, 750 & 1000 ml/minute
- Supplied with flexible tube for connection to your application, 9V battery and Certificate of Calibration

<b>Flow (ml/min)</b>	<b>Display Resolution (ml/min)</b>
0.50-9.99	0.01
10.0	99.9
100-1000	1.0